INSTRUCTOR'S MANUAL For Use With PRACTICAL RADIOGRAPHIC IMAGING Eighth Edition

Eighth Edition

INSTRUCTOR'S MANUAL For Use With PRACTICAL RADIOGRAPHIC IMAGING

By **QUINN B. CARROLL, M.Ed., R.T.**

Director Radiography Department Midland College Midland, Texas



CHARLES C THOMAS • PUBLISHER, LTD. Springfield • Illinois • U.S.A.

Published and Distributed Throughout the World by

CHARLES C THOMAS • PUBLISHER, LTD. 2600 South First Street Springfield, Illinois 62704

This book is protected by copyright. No part of it may be reproduced in any manner without written permission from the publisher. All rights reserved.

$\ensuremath{\mathbb{C}}$ 2007 by CHARLES C THOMAS • PUBLISHER, LTD.

ISBN 978-0-398-07710-x (paper)

With THOMAS BOOKS careful attention is given to all details of manufacturing and design. It is the Publisher's desire to present books that are satisfactory as to their physical qualities and artistic possibilities and appropriate for their particular use. THOMAS BOOKS will be true to those laws of quality that assure a good name and good will.

Printed in the United States of America R-3

CONTENTS

Pag
ADDITIONAL ASSIGNMENT QUESTIONS: CALCULATIONS 3
Chapter 5
Chapter 6
Chapter 9
Chapters 10 & 11
Chapter 13
Chapter 14
Chapter 18
Chapter 20
Chapter 26
Chapter 27
Chapter 31
ANSWERS TO CALCULATIONS ASSIGNMENT
QUESTIONS
ANSWERS TO CHAPTER REVIEW QUESTIONS19
MULTIPLE CHOICE QUESTIONS BANKS FOR
EACH CHAPTER47
ANSWERS TO MULTIPLE CHOICE QUESTIONS161
SUGGESTED LABORATORY EXERCISES FOR
EACH CHAPTER169
General Instructions for Labs
Lab Practice #1: Determining Density and Contrast Changes 174
Lab #1: mAs
Lab #2: Adequate Penetration
Lab #3: Kilovoltage-Peak
Lab #4: 15 Percent Rule for kVp180
Lab #5: Machine Phase
Lab #6: Field Size Limitation
Lab #7: Grids
Lab #8: Grid Cut-Off
Lab Practice #2: Measuring Sharpness of Recorded Detail,
Magnification and Distortion
Lab #9: Intensifying Screen Emission
Lab #10: Intensifying Screens

Lab #11:	Focal Spot Size
Lab #12:	Anode Heel Effect196
Lab #13:	Focus-Film Distance
Lab #14:	Object-Film Distance
Lab #15:	SID/OID Ratio
Lab #16:	Beam-Part-Film Alignment
Lab #17:	Timer Quality Control Tests
Lab #18:	Milliamperage Quality Control Tests205
Lab #19:	Kilovolt-Peak and Half Value Layer Quality
Contro	l Tests
Lab #20:	Focal Spot Quality Control Tests
Lab #21:	Beam Alignment and Field Size Tests
Lab #22:	Sensitometry
Lab #23:	Processor Control Charts
Lab #24:	Personal Repeat Analysis

INSTRUCTOR'S MANUAL For Use With PRACTICAL RADIOGRAPHIC IMAGING Eighth Edition

ADDITIONAL ASSIGNMENT QUESTIONS: CALCULATIONS

Chapter 5-mAs Calculations Bank (oral or written)

-	
I. $mA \times s = mAs $ (5 seconds each)	II. $mAs = mA \times s $ (8 seconds each)
	MUST BE FRACTIONS, NOT DECIMALS:
1. $50 \times 1/20 =$	16. $10 \text{ mAs} = 50 \text{ mA x}$
2. $200 \times 1/5 =$	17. 3.3 mAs – 100 mA x
3. $100 \times 1/12 =$	18. $1.6 \text{ mAs} = 100 \text{ mA x}$
4. $200 \times 1/8 =$	19. $15 \text{ mAs} = 300 \text{ mA x}$
5. $200 \times 7/20 =$	20. $45 \text{ mAs} = 300 \text{ mA x}$
6. $100 \times 3/20 =$	21. $8.3 \text{ mAs} = 100 \text{ mA x}$
7. $300 \times 3/20 =$	22. $60 \text{ mAs} = 300 \text{ mA x}$
8. $300 \times 1/8 =$	23. 80 mAs = 400 mA x
9. $100 \times .05 =$	24. $240 \text{ mAs} = 400 \text{ mA x}$
10. $200 \times .15 =$	25. $66 \text{ mAs} = 100 \text{ mA x}$
11. $300 \times .25 =$	26. 4 mAs = 100 mA x
12. $200 \times .07 =$	27. 5 mAs = 150 mA x
13. $400 \times .02 =$	28. 180 mAs = 300 mA x
14. $200 \times .02$	29. 12.5 mAs = 100 mA x
15. $50 \times .4 =$	30. 160 mAs = 200 mA x
greater than the original?	
Fill in the kVp in the following sity with the original technique,	problems which would maintain equal denlisted first:
11. 200 mA 1/8 sec. 60 kVp	13. 50 mA 1/60 sec. 50 kVp
100 mA 1/2 sec kVp	400 mA 1/120 sec kVp
12. 400 mA 1/20 sec. 90 kVp	14. 300 mA 1/30 sec. 70 kVp
300 mA 1/60 sec kVp	150 mA 1/10 sec kVp